

Chromatography Supplies




Vials and Septa-Caps Kit • Long Life D2 Lamps • Hollow Cathode Lamps




	Part No.	Item Name	Pack Qty.
		9mm Screw Top Vials	
A	ST-09-2067	1.5ml Clear Glass Screw type vials, w/o write-on	100
B	ST-09-2068	1.5ml Clear Glass Screw type vials, with write-on	100
C	ST-09-2070	1.5ml Amber Glass Screw type vials, w/o write-on	100
D	ST-09-2069	1.5ml Amber Glass Screw type vials, with write-on	100
E	ST-09-0932	0.3 ml (300 µl) PP micro Vial	100
F	ST-PP-1205	1.5ml PP Screw type vials	100
G	ST-09-0620	1.1ml Clear Glass HIGH RECOVERY Screw type vials, w/o write-on	100



	Part No.	Item Name	Pack Qty.
		9mm Cap and Septa	
A	SCSPB-09-0730	PP Blue Colour Open Top Screw Cap, 6mm hole with Red PTFE / White Silicon BONDED septa for 9-425 vials	100
B	SCSPB-09-2074	PP Blue Colour Open Top Screw Cap, 6mm hole with Red PTFE / White Silicon SLIT BONDED septa for 9-425 vials	100
C	SCSPB-09-1534	PP Blue Colour Open Top Screw Cap, 6mm hole with White PTFE / Beige Silicon SLIT BONDED septa for 9-425 vials	100
D	SCSP-09-0730	PP Blue Colour Open Top Screw Cap, 6mm hole with Red PTFE / White Silicon septa for 9-425 vials	100
E	SCSP-09-2074	PP Blue Colour Open Top Screw Cap, 6mm hole with Red PTFE / White Silicon with SLIT septa for 9-425 vials	100



	Part No.	Item Name	Pack Qty.
		Inserts	
A	MI-09-0357	0.1ml, 31 x 6mm micro insert (For 9mm vial)	100
B	MI-09-0865	0.1ml, 31 x 6mm micro insert with assembled plastic spring (For 9mm vial)	100
C	MI-09-0866	0.2ml, 31 x 6mm flat bottom micro insert (For 9mm vial)	100



	Part No.	Item Name	Pack Qty.
		11mm Crimp Top Vials, Cap and Septa	
A	CT-11-4491	1.5ml Clear Glass Crimp type vials, w/o write-on	100
B	CT-11-4492	1.5ml Clear Glass Crimp type vials, with write-on	100
C	CT-11-4493	1.5ml Amber Glass Crimp type vials, with write-on	100
D	CT-11-4494	1.5ml Amber Glass Crimp type vials, w/o write-on	100
E	CC-11-4500	11mm Silver Aluminium Crimp Caps with Center hole, with White Silicon/Red PTEF Septa	100



	Part No.	Item Name	Pack Qty.
		11mm Snap Top Vials, Cap and Septa	
A	ST-11-1232	1.5ml Clear Glass Snap Top vials, w/o write-on	100
B	ST-11-1232W	1.5ml Clear Glass Snap Top vials, with write-on	100
C	ST-11-1232A	1.5ml Amber Glass Snap Top vials, w/o write-on	100
D	ST-11-1232WA	1.5ml Amber Glass Snap Top vials, with write-on	100
E	SCSP-11-4500	11mm Clear Snap Cap, 6mm hole with Red PTFE septa/ White Silicon Septa	100



	Part No.	Item Name	Pack Qty.
		8mm Screw Top Vials, Cap and Septa	
A	ST-08-4427	1.5ml Clear Glass Screw type vials, w/o write-on	100
B	ST-08-4428	1.5ml Clear Glass Screw type vials, with write-on	100
C	ST-08-4429	1.5ml Amber Glass Screw type vials, w/o write-on	100
D	ST-08-4430	1.5ml Amber Glass Screw type vials, with write-on	100
E	SCSP-08-4437	PP Black Colour Screw Cap, with Red PTFE /White Silicon septa	100



	Part No.	Item Name	Pack Qty.
		13mm Screw Top Vials, Cap and Septa	
A	ST-13-0222	4ml Clear Glass Screw Top Vial, w/o write-on	100
B	ST-13-1335	4ml Clear Glass Screw Top Vial, with write-on	100
C	ST-13-0280	4ml Amber Glass Screw Top Vial, w/o write-on	100
D	ST-13-1336	4ml Amber Glass Screw Top Vial, with write-on	100
E	SCSP-13-4437	PP Black Colour Open Top Screw Cap, 8.5mm hole with Red PTFE/ White Silicone septa	100

	Part No.	Item Name	Pack Qty.
		1ml Shell vial	
A	SV-01-0513	1ml Clear Glass Shell Vial with PE Plug Cap	100
B	SV-PP-2194	1ml PP Shell Vial with PE Plug Cap	100



	Part No.	Item Name	Pack Qty.
		20ml Headspace Vials, Cap and Septa	
A	PE-0C-0297	20ml Clear Glass Flat Bottom Vial for PE	100
B	PE-0C-0297 (R)	20ml Clear Glass Round Bottom Vial for PE	100
C	UN-RB-0873	20ml Clear Glass Round Bottom Vial	100
D	UN-FB-0796	20ml Clear Glass Flat Bottom Vial	100
E	ALNS-STAG-6579	20mm, Silver Aluminium Crimp Cap with Silicon White /PTFE Beige Dark (Tan) Septa	100
F	MG-08-6579	20mm, Gold Magnetic Crimp Cap with Silicon White /PTFE Beige Dark (Tan) Septa (caps & Septa packed separate)	100
G	BM-08-6579	20mm, Bimetallic Crimp Cap with Silicon White /PTFE Beige Dark (Tan) Septa (caps & Septa packed separate)	100
H	ALNS-STAG-0638	20mm, Silver Aluminium Crimp Cap, Silicon White /PTFE Beige (Tan) Septa	100
I	MG-08-0638	20mm, Gold Magnetic Crimp Cap with Silicon White /PTFE Beige (Tan) Septa (caps & Septa packed separate)	100

	Part No.	Item Name	Pack Qty.
		20ml Head Space screw type vial, Cap and Septa	
A	ST-20-1307	20 ml Clear Glass Round Bottom Screw type Vial	100
B	ST-20-1311	20 ml Amber Glass Round Bottom Screw type Vial	100
C	SCSP-20-1309	18 mm Silver Magnetic Screw Cap with Silicon White /PTFE Blue Septa	100

	Part No.	Item Name	Pack Qty.
		EPA KIT (Certified)	
A	KT-40-0402	40ml Clear Glass Screw Neck Vial with PP White Colour Open Top Screw Cap with Silicon natural /PTFE beige BONDED Septa	100

	Part No.	Item Name	Pack Qty.
		Cap and Septa	
A	SCSP-24-0842	PP White Colour Open Top Screw Cap with Silicon natural /PTFE beige BONDED Septa (EPA -Quality)	100

	Part No.	Item Name	Pack Qty.
		Crimper and Decapper	
A	CRIMP-1	Manual Crimper 11mm	1
B	DECAP-11	Manual Decapper 11mm	1
C	CRIMP-20	Manual Crimper 20mm	1
D	DECAP-20	Manual Decapper 20mm	1



Guaranteed quality long life hplc lamps

Deuterium 2000 Hour Lamps
for HPLC UV Detection

Key Points

- OEM equivalent lamps
- Same bare lamps as used to produce
- OEM lamp
- Guarantees and warranties match those offered by OEM

Agilent

1100 1200 1260 1290

Cyberlab Part Number	Description
D2-AGI-VWD	D2 Lamp for Agilent 1100/1200 VWD/UV equivalent to G1314-60100/G1314-60101
D2-AGI-DAD	D2 Lamp for Agilent 1100/1200 DAD/PDA equivalent to P/N 5181-1530/2140-0813/2140-0820
D2-AGI-1260	D2 Lamp for Agilent 1260/1290 equivalent to 5190-0917

Waters

996 2996 2487 Alliance

Cyberlab Part Number	Description
D2-WAT-2487	D2 Lamp for Waters Alliance 2487 UV/VWD equivalent to P/N WAS081142
D2-WAT-2996	D2 Lamp for Waters Alliance 2996 PDA/DAD equivalent to P/N WAT052586/WAT057760
D2-WAT-UPLC	D2 Lamp for Waters ACQUITY UHPLC 2489/2998 equivalent to P/N 201000186

Shimadzu

SPD10A AVP AVVP M10AVP 20A 20AV 2010

Cyberlab Part Number	Description
D2-SHI-2010	D2 Lamp for Shimadzu HPLC LC-2010 equivalent to P/N 228-37401-00
D2-SHI-SPD	D2 Lamp for Shimadzu SPD10A/VP equivalent to P/N 228-34016-02
D2-SHI-UVAA	D2 Lamp for Shimadzu AA/UV Spectrophotometers equivalent to P/N 062-65055-05

Shimadzu Prominence Lamp

Cyberlab Part Number	Description
D2-SHI-NEXPROD	D2 Lamp for Shimadzu NEXERA (LC 2040/2040C) & Prominence I (LC2030/2030C) equivalent to P/N 228-55626-01



Hollow Cathode Lamps

The Cyberlab range of high performance Single Element Hollow Cathode lamps are available in 37mm (1.5") or 50mm (2") in both Coded and Uncoded version.

Hollow Cathode Lamps for AA Spectrophotometers

- Hollow Cathode Lamps 37mm (1.5") Diameter Single Element
- Hollow Cathode Lamps 37mm (1.5") Diameter Multiple Element's
- Hollow Cathode Lamps THERMO FISHER / UNICAM 37mm (1.5") Diameter Single Element 7Pin Base
- Hollow Cathode Lamps AGILENT / VARIAN 37mm (1.5") Diameter Single Element 4Pin Base
- Hollow Cathode Lamps 51mm (2") Diameter Single Element to suit Perkin Elmer
- Hollow Cathode Lamps 51mm (2") Diameter Single Element to suit Perkin Elmer ANALYST
- Hollow Cathode Lamps 51mm (2") Diameter Multiple Element's to suit Perkin Elmer

List of Elements

Element	Symbol	Element	Symbol	Element	Symbol	Element	Symbol
Aluminium	Al	Gadolinium	Gd	Nickel	Ni	Strontium	Sr
Antimony	Sb	Gallium	Ga	Niobium	Nb	Tantalum	Ta
Arsenic	As	Germanium	Ge	Osmium	Os	Tellurium	Te
Barium	Ba	Gold	Au	Palladium	Pd	Terbium	Tb
Beryllium	Be	Hafnium	Hf	Phosphorus	P	Thallium	Tl
Bismuth	Bi	Holmium	Ho	Platinum	Pt	Thorium	Th
Boron	B	Indium	In	Potassium	K	Thulium	Tm
Cadmium	Cd	Iridium	Ir	Praseodymium	Pr	Tin	Sn
Calcium	Ca	Iron	Fe	Rhenium	Re	Titanium	Ti
Cerium	Ce	Lanthanum	La	Rhodium	Rh	Tungsten	W
Caesium	Cs	Lead	Pb	Rubidium	Rb	Uranium	U
Chromium	Cr	Lithium	Li	Ruthenium	Ru	Vanadium	V
Cobalt	Co	Lutetium	Lu	Samarium	Sm	Ytterbium	Yb
Copper	Cu	Magnesium	Mg	Scandium	Sc	Yttrium	Y
Deuterium	D	Manganese	Mn	Selenium	Se	Zinc	Zn
Dysprosium	Dy	Mercury	Hg	Silicon	Si	Zirconium	Zr
Erbium	Er	Molybdenum	Mo	Silver	Ag		
Europium	Eu	Neodymium	Nd	Sodium	Na		



Varian Coded & Unicam Coded Lamps are also available