

# PRODUCTS



## FILTRATION CUPS:

- Pre-sterilized, single-use devices for efficient microbial analysis in pharmaceutical API and formulation industries. Individually packed for safety and convenience.
- User-friendly design: 100 ml detachable funnel and sealable base for easy sample collection and lab transport. Boost efficiency with our filtration cups.

## SPECIMEN CONTAINER:

- Available in both sterile and non-sterile options, suitable for sample collection.
- Perfect for collecting, storing, and transporting samples with ease.
- Sterile containers undergo gamma irradiation for safety.
- Container made of polystyrene, cap crafted from high-density polyethylene.

## DISPOBAGS:

- Suitable for sterilizing glassware or plastic wares.
- Autoclavable bags designed for waste disposal, ideal for disposing of pathological, clinical, or contaminated materials.

## SPOON/SPATULA:

- Available in both sterile and non-sterile options.
- Individually packed, featuring a flat surface and small scoop.
- Ideal for weighing dry powders.
- Sterile variants undergo gamma irradiation for safety.

## Type 1/DM Water:

- **Ultra-Pure Water** : High purity and UV transmission for column and detector protection, ensuring reproducibility.
- **HPLC Water**: Ideal for reproducible separations due to high UV transmittance, low particle count, acidity, alkalinity, and evaporation residue.
- **DM Water**: Suitable for various lab applications, commonly used in regulated and demanding settings.



## Address:

Kamgo Towers, Kamgo Technologies India Pvt. Ltd., 141, 9th Link Road, Nehru Nagar, Industrial Estate, Chennai 600041